



Protocol NT001:

Neuronal Tract Tracing in Fixed Tissue Using NeuroVue[®] Dye-Coated Filters

Overview of Labeling Strategy (see page 2 for more detailed instructions)

- Fix tissue in 4% buffered formaldehyde.
- Initiate labeling by inserting NeuroVue micro-strip(s) into nerve tract(s) to be traced. The highly lipophilic NeuroVue dyes transfer from the micro-strip into nerve cell membranes and diffuse along the lipid bilayer in both directions from the insertion site (anterograde and retrograde labeling).
- Incubate tissue in 4% phosphate buffered formaldehyde at 37°C.
- Monitor the progress of dye diffusion using light microscopy and/or fluorescence microscopy.
- When dye(s) have reached the region(s) to be studied, remove NeuroVue micro-strip(s) and prepare whole mounts or tissue sections for fluorescence imaging.

Advantages of NeuroVue technology

- Different fibers can be traced in the same specimen by using fluorescent NeuroVue dyes that excite and emit in the green, orange, red and/or far red (<http://www.mtarget.com/product/neurovue.aspx>).
- Neuronal connections can be studied in embryos lacking receptors needed for neuronal identification as well as in juveniles and adults (*Gurung & Fritsch, J Comp Neurol 479:309-327, 2004; Morris et al., Brain Res 1091:186-199, 2006; Hsieh & Cramer, J Comp Neurol 497:589-599, 2006*).
- Use of dye-coated filters allows more precise positioning than is possible with crystals or oils, avoids tissue damage caused by high pressure microinjection, and provides sharp high resolution images of both afferent and efferent fibers arising at the point of filter insertion (*Fritsch et al., Brain Res Bull 66:249-258, 2005*).
- Use of NeuroVue dyes reduces the complexity of labeling procedures because these dyes have been selected to have similar diffusion rates, allowing simultaneous or near-simultaneous application of different colors in most cases (*Fritsch et al., Brain Res Bul, 66:249-258, 2005*).

Equipment and Supplies Needed

Item	Micro-strip Preparation	Specimen Fixation	Labeling/ Diffusion Monitoring	Tissue Mounting
NeuroVue Filter square (1cm x 1cm)	X			
Protective gloves	X	X	X	X
Microscissors*	X			
Forceps	X	X	X	X
95% ethyl alcohol	X	X		
Dissecting microscope	X		X	
4% formaldehyde [#] ; pH 7.4 in phosphate buffered saline (PBS)		X	X	
Refrigerator or cold room (4 °C.)		X		
NeuroVue Micro-strips (see Step 1, page 2)			X	
Screw-cap vials for specimen storage/incubation		X	X	
37 °C incubator or oven			X	
15% gelatin in PBS (or other non-paraffin embedding material)				X
Glycerol (or other mounting medium free of organic solvents)				X
Microscope slides and coverslips				X
Fluorescence microscope with confocal and/or epi-illumination)				

* World Precision Instruments #500086, (<http://store.wpiinc.com>) or equivalent

[#] Prepared from ultrapure formaldehyde solution (e.g., Polysciences, #04018, www.polysciences.com/shop/ or equivalent) or reagent grade paraformaldehyde powder (e.g., Sigma P-6148, www.sigmaaldrich.com or equivalent).

PROCEDURES

The examples below, kindly provided by Drs. Bernd Fritzsich and Lucy Feng (Creighton University), illustrate applications of NeuroVue® dyes for tracing various sensory nerves of the head and ear to and from the brain. For additional examples in a variety of developmental systems, see the NeuroVue Bibliography below.

1. Preparation of NeuroVue Filter Micro-strips

- a) Put on a clean pair of protective gloves.
- b) Agitate scissors and forceps in alcohol to remove any dye residues. Air dry completely before use.

Note: To avoid cross-contamination among dyes, it is preferable to use a separate set of instruments for each dye being handled in Steps 1, 3 and 4.

Note: It is important to avoid residual alcohol on tools used to handle the NeuroVue filters or micro-strips, especially in Steps 1 and 3. Contact with alcohol (or other organic solvents) can potentially alter the standardized dye concentration coated on the filter. This in turn may cause increased variability in diffusion times from sample to sample or study to study.

- c) Using a dissecting microscope, hold the NeuroVue dye coated filter with forceps and use microscissors to cut into triangular pieces approximately 0.3 -1 mm on a side (**Figure 1**; to view details, enlarge image electronically).

Note: Shape and size of pieces can be adjusted depending upon the size of the tissue specimen to be labeled and the number of colors to be used simultaneously.

- d) After cutting micro-strips, agitate scissors and forceps in alcohol to remove any dye residues, then air dry completely.




Figure 1. Preparation of NeuroVue Red micro-strips for use in tissue labeling. Using a dissecting microscope, microscissors (at left) are used to cut small triangles from a 1 cm x 1 cm coated filter square (Product #FS-1002). Scale bar = 2 mm (~25X magnification).

2. Preparation of Tissue Specimens for Labeling with NeuroVue Micro-strips

- a) Put on a pair of protective gloves
- b) Perfusion fix embryos with 4-10% formaldehyde in 0.1M phosphate buffer (pH 7.4).

Note: In order to minimize background tissue fluorescence caused by autofluorescent impurities in the fixative it is important that fixative be prepared using ultrapure methanol free formaldehyde (e.g. Polysciences catalog #04018) or paraformaldehyde.

- c) Store fixed specimen at 4°C in 4% neutral buffered formaldehyde for at least 3 days (or up to 1 year)

Note: Prolonged fixation (>1 month) results in greater tissue crosslinking and slower dye diffusion.

- d) Rinse forceps to be used for manipulation with alcohol to remove any dye residues.
- e) For brain, remove and pin down the head. If needed, dissect to expose application site (see *Fritzsich Jackson Lab Presentation 2005* for examples of “inside-out” labeling of the inner ear). For more rigid tissues, make incision(s) as needed to allow insertion of filter micro-strip(s) and note a landmark to find the cut later.

Note: To avoid spread of dye through unwanted neuronal profiles through peripheral anastomoses, cut nerve(s) for which labeling is not desired. For example, cut the facial nerve root to avoid cross labeling of the trigeminal nerve, and vice versa.

3. Insertion of NeuroVue Micro-strips in Tissue to Initiate Labeling

- a) Put on a pair of protective gloves
- b) Agitate forceps in alcohol before and after insertion of each micro-strip. Air dry completely before use.

Note: To avoid cross-contamination, use a separate forceps for each different NeuroVue dye.

Note: It is important to avoid residual alcohol on tools used to handle the NeuroVue filters or micro-strips, especially in Steps 1 and 3. Contact with alcohol (or other organic solvents) can potentially alter the standardized dye concentration coated on the filter. This in turn may cause increased variability in diffusion times from sample to sample or study to study.

- c) Pick up NeuroVue micro-strip with forceps and insert in tissue site.
- d) Insert additional NeuroVue micro-strips as needed (**Figure 2**).

Note: In soft tissue, such as brain, the filter can be directly inserted with the sharp (tip) end first.

Note: In more rigid tissues, re-locate the landmark(s) noted during tissue preparation (Step 3) and insert the micro-strip, using its tip to assist in finding the incision.

Note: Do NOT use a dissecting needle to push filter into tissue. This can cause inaccurate placement and also inadvertent labeling of non-target fibers.

- e) Place tissue with inserted micro-strip(s) in a closed jar with 4% buffered formaldehyde and incubate at 22°C – 37°C, monitoring periodically as described in Step 4 below, until neuronal profiles of interest are filled.

Note: Increasing incubation temperature decreases structural stability of fixed tissue but decreases diffusion times required for filling of neuronal profiles.

4. Monitoring Dye Diffusion

- a) Estimation of diffusion times needed

Diffusion times must generally be determined empirically for each laboratory's neuronal profile(s) of interest since they depend on a number of variables:

- i) Distance over which neuronal profile(s) are to be traced (**Table 1**);
- ii) Dye concentration (standardized for each NeuroVue dye);
- iii) Incubation temperature (increasing temperature decreases diffusion time but may reduce structural preservation; temperatures substantially above 37°C may reduce dye stability);
- iv) Extent of tissue cross linking (higher fixative concentration and/or longer times in fixative will increase diffusion times).

- b) Checking Diffusion Distance

- i) Put on a pair of protective gloves
- ii) Use a dissecting microscope to observe how far the dye front has traveled from the point of micro-strip insertion. All of the NeuroVue dyes absorb strongly in the visible: NeuroVue Green/ Emerald/Jade appear yellow/orange, NeuroVue Red appears magenta/red and NeuroVue Maroon appears cyan/blue (**Figure 3**).

Note: Fluorescence detection of NeuroVue dyes is much more sensitive than detection using absorbance. Thus, use of absorbance for monitoring (Figure 3) will underestimate the true extent of neuronal profile filling. If dye front extends to the final target region when observed using a dissecting microscope, diffusion time may have been too long.




Figure 2. Placement of NeuroVue micro-strips for multicolor neurotracing. Lateral view of murine head (embryonic day 12.5); with micro-strips placed to obtain central projections of NeuroVue Red labeled trigeminal nerve, NeuroVue Maroon labeled facial nerve and NeuroVue Green labeled glossopharyngeal nerve. The eye is visible as a brown spot at left (anterior). Same magnification as in Figure 1 (~25X).

Table 1. Approximate Diffusion Times for Filling of Murine Ear Projections to and from the Brain

Developmental Stage	Time @ 37°C.
Embryonic day 11.5 (E11.5)	24 – 36h
Embryonic day 13.5 (E13.5)	48 – 60h
Embryonic day 16.5 (E16.5)	72 – 80h
Newborn (P0)	80 – 96h





Figure 3. Monitoring diffusion distance using NeuroVue dye absorbance. After incubation for 36 h at 37°C, diffusion in all directions from the point of micro-strip insertion is readily visualized using a dissecting microscope (same specimen as in Figure 2; magnification ~25X).

Note: Dye concentrations high enough to give visible color may also cause absorption quenching of emitted fluorescence.

Note: Like other membrane dyes used for neurotracing, the NeuroVue dyes will also label muscle fibers or any other macroscopic structures surrounded by continuous lipid bilayers (Figure 4). Care must therefore be taken in Step 3 to insert the micro-strip in a location that maximizes labeling of target neurons and minimizes labeling of non-target structures.

Figure 4. Monitoring proper placement of NeuroVue micro-strips. Correct insertion in murine spinal cord labels only intercostal nerves after incubation for 10 days at 37°C (not shown). Suboptimal insertion can also label intercostal muscle fibers to varying extents (here, strong labeling for NeuroVue Orange, moderate for NeuroVue Maroon). Scale bar = 2 mm.




5. Preparation of Labeled Tissue For Imaging

Note: For the following procedures, wear protective gloves. Agitate tools used for tissue manipulation in alcohol and air dry between steps to avoid dye cross contamination during the dissection process.

a) Whole mounts or surface preparations

Figure 5 below illustrates the major steps involved in preparing and imaging whole brain mounts for the type of specimen shown in Figure 2. Use of such preparations in combination with confocal microscopy is the preferred method for imaging patterns of innervation in brain or peripheral structures such as ear, cochlea, etc. (see **Figure 6** below)



b) Serial sections

Use of such preparations is essential when studying larger specimens (e.g., older stages of brain development) and provides additional information for complex peripheral structures such as the ear. Embedding in 15% gelatin, post-fixation for at least 48 h in 4-10% buffered formaldehyde, preparation

of serial vibratome sections and mounting in glycerol (*Gurung & Fritsch, J Comp Neurol 479:309-327, 2004*) is the preferred method of tissue preparation when whole mounts or surface mounts do not allow adequate visualization of tissue regions of interest.

Note: The following precautions have been found to minimize transcellular dye spreading and maximize ability to maintain good resolution of labeled fibers when preparing serial sections from gelatin-embedded tissue:

- Prepare sections in chilled (4°C) buffer and keep them in chilled buffer until just before mounting.
- Image quality may be reduced if sections are mounted in glycerol and held at room temperature for longer than 1 hour before viewing. Therefore delays between mounting and imaging should be minimized by mounting and coverslipping only a limited number of sections (8-12) in 100% glycerol at any one time, and viewing or imaging them as soon as practical.
- If sections must be held in glycerol for longer than 1 hour, store at 4°C. Some image degradation will be seen in cut profiles at the tissue surface but overall image quality will be good even after overnight storage at 4°C.

Note: Frozen sections are less optimal due to greater dye leakage during the sectioning process (*von Bartheld et al., J. Histochem Cytochem 38:725-733, 1990; Köbbert et al., Progr Neurobiol 62: 327-351, 2000*).

Note: As with other membrane dyes used for neurotracing, paraffin embedding or use of organic solvents should be avoided when using the NeuroVue dyes. Such solvents extract membrane lipids and associated dyes, reducing or eliminating ability to detect labeled nerve fibers.

6. Imaging of Whole Mounts or Tissue Sections

As of January 2008, NeuroVue dyes are available that enable 3 – 5 color neuronal profiling in combination with commonly used genetic markers. **Table 2** briefly summarizes the spectral properties of the commercially available NeuroVue dyes. Full excitation and emission spectra, appropriate laser lines and filter combinations for confocal imaging, and suggested excitation-emission filter sets for epifluorescence imaging of each of these dye may be found in the individual product data sheets.

Which dyes are optimal depends on the instrument configuration available and the length of time required for diffusion from point of filter insertion to the region of interest (**Table 3**). For some studies with tissue which exhibits high levels of green/yellow autofluorescence background, it is recommended that the green channel be reserved for genetic markers and the longer wavelength channels for the NeuroVue dyes.

Note: Due their very long red fluorescence emissions, most people cannot see NeuroVue Maroon or NeuroVue Burgundy emissions by eye. Detection by camera will be more sensitive than with the unaided eye.

Product Number	Product Name	Excitation maximum	Emission Maximum	Ext. coeff. (M ⁻¹ cm ⁻¹)
FS-1001	NeuroVue Maroon	647 nm	667 nm	222,368
FS-1002	NeuroVue Red	567 nm	588 nm	120,000
FS-1003	NeuroVue Orange	550 nm	570 nm	139,800
FS-1005	NeuroVue Burgundy	683 nm	707 nm	192,832
FS-1006	NeuroVue Jade	478 nm	508 nm	87,850

Type of Study	Dye Combination	Instrument requirement
3 color neurotracing Medium term (≤ 5 days)	NeuroVue Jade NeuroVueRed, NeuroVue Maroon	For standard epifluorescence microscopes and confocal systems having 488nm, 568nm and 633/647nm excitation
3 color neurotracing Medium term (≤ 5 days)	NeuroVue Jade, NeuroVue Orange NeuroVue Maroon	For standard systems with 543nm excitation instead of 568nm
4 color neurotracing Long term studies (3-4 wk) ¹	NeuroVue Orange, NeuroVue Red NeuroVue Maroon NeuroVue Burgundy	For systems with spectral detection.
5 color neurotracing Medium term (≤ 5 days)	NeuroVue Jade NeuroVue Orange, NeuroVue Red NeuroVue Maroon NeuroVue Burgundy	Color unmixing required for Orange/Red and Maroon/Burgundy combinations

¹ Compatible with green fluorescent genetic markers such as eGFP or photoactivated X-Gal reaction product BCI [*Matei, Brain Res Bull 70:33-43, 2006*]

Note: New NeuroVue family members are also in development at MTTI. If you have special needs that are not met by the products listed in Table 2, please contact us via phone, fax or e-mail to inquire.

As shown in the examples below, the NeuroVue dyes can be used for neuroanatomical tracing at many levels of resolution, in combination with one another or with genetic markers (e.g., eGFP) specific to neuronal subtypes of interest.

Figure 6 illustrates confocal imaging of facial and cranial nerve tracts using a whole brain mount.




Figure 7 illustrates confocal imaging of innervation in the murine inner ear using a surface preparation of the cochlea. Note that even very thin Type II inner hair cells (IHC), which are normally quite difficult to visualize, are clearly visible at the single fiber level (right panel).




Figure 8 illustrates confocal imaging of the interactions between radial fibers in the murine inner ear and supporting cells in the organ of Corti, again using a surface preparation of the cochlea.





Figure 9 illustrates the use of two color neurotracing to optically dissect anatomically overlapping fibers.



For further examples of how the NeuroVue dyes have been used to trace neuronal interactions in a wide variety of systems, see the NeuroVue Bibliography below, watch for emerging publications on our website (www.mtarget.com), and send us yours to add to the list!

REFERENCES CITED

- Fritzsich B, Muirhead KA, Feng F, Gray BD and Ohlsson-Wilhelm BM. 2005. Diffusion and imaging properties of three new lipophilic tracers, NeuroVue Maroon™, NeuroVue Red™ and NeuroVue Green™ and their use for double and triple labeling of neuronal profile. *Brain Res Bull* 66: 249-258.
- Gurung B, Fritzsich B. 2004. Time course of embryonic midbrain and thalamic auditory connection development in mice as revealed by carbocyanine dye tracing. *J Comp Neurol* 479:309-327.
- Hsieh CY, Cramer KS. 2006. Deafferentation Induces Novel Axonal Projections in the Auditory Brainstem After Hearing Onset. *J Comp Neurol* 497: 589-599.
- Köbber C, Apps R, Bechmann I, Lanciego JL, Mey J, Thanos S. 2000. Current concepts of neuroanatomical tracing. *Progress in Neurobiology* 62: 327-351.
- von Bartheld CS, Cunningham DE, Rubel EW. 1990. Neuronal tracing with Dil: decalcification, cryosectioning, and photoconversion for light and electron microscopic analysis. *J Histochem Cytochem* 38:725-733.

NeuroVue® BIBLIOGRAPHY

- Beisel, K.W., Wang-Lundberg, Y., Maklad, A. and Fritzsich, B. 2005. Development and evolution of the vestibular sensory apparatus of the mammalian ear. *J Vestibular Research* 15:225-241. **NeuroVue Maroon** (previously PTIR271; see Hellard, *Dev. Biol.* 2004, below, for methods).
- de Caprona MD, Beisel KW, Nichols DH, Fritzsich B. 2004. Partial behavioral compensation is revealed in balance tasked mutant mice lacking otoconia. *Brain Res Bull* 64:289-301. **Both NeuroVue Maroon** (previously PTIR271) and **NeuroVue Red** (previously PTIR278) were used in Figure 8 (B. Fritzsich, personal communication).
- Fritzsich B. 2003. Development of inner ear afferent connections: forming primary neurons and connecting them to the developing sensory epithelia. *Brain Res Bull* 60:423-433. **NeuroVue Maroon** (previously PTIR271) and **NeuroVue Green** (previously PTIR281) were used in Figure 4a, 4c, and 4d; **NeuroVue Maroon** and PKH26 were used for cover image (B. Fritzsich, personal communication)
- Fritzsich B, Tessarollo L, Coppola E, Reichardt LF. 2004. Neurotrophins in the ear: their roles in sensory neuron survival and fiber guidance. *Prog Brain Res* 146:265-278. **NeuroVue Maroon** (formerly PTIR271) and *Dil* were used in Figure 2 (B. Fritzsich, personal communication)
- Fritzsich B, Muirhead KA, Feng F, Gray BD, Ohlsson-Wilhelm BM. 2005. Diffusion and imaging properties of three new lipophilic tracers, NeuroVue Maroon, NeuroVue Red and NeuroVue Green and their use for double and triple labeling of neuronal profile. *Brain Res Bull* 66:249-258. **NeuroVue Maroon, NeuroVue Red, NeuroVue Green**
- Fritzsich B, Matei VA, Nichols DH, Bermingham N, Jones K, Beisel KW, Wang VY. 2005. Atoh1 null mutants show directed afferent fiber growth to undifferentiated ear sensory epithelia followed by incomplete fiber retention. *Dev Dyn*, 233: 570-583. **NeuroVue Maroon** (previously PTIR271), **NeuroVue Red** (previously PTIR278)
- Fritzsich B, Pauley S, Matei V, Katz DM, Xiang M, Tessarollo L. 2005. Mutant mice reveal the molecular and cellular basis for specific sensory connections to inner ear epithelia and primary nuclei of the brain. *Hear Res*, 206: 52-63. **NeuroVue Maroon** (previously PTIR271); see Hellard, *Dev Biol* 2004, below, for methods)
- Fritzsich B, Jackson Lab Presentation, 2005:
http://www.biomedsci.creighton.edu/facilities/nccb/media/Jackson_lab_presentation.ppt
NeuroVue Green (previously PTIR281); **NeuroVue Red** (previously PTIR278); **NeuroVue Maroon** (previously PTIR271)

- Gu C, Rodriguez ER, Reimert DV, Shu T, Fritsch B, Richards LJ, Kolodkin AL, Ginty DD. 2003. Neuropilin-1 Conveys Semaphorin and VEGF Signaling during Neural and Cardiovascular Development. *Dev Cell* 5:45-57. **NeuroVue Maroon** (formerly PTIR271) rather than PKH26 was used in Figure 3 (B. Fritsch, personal communication)
- Gurung B, Fritsch B. 2004. Time course of embryonic midbrain and thalamic auditory connection development in mice as revealed by carbocyanine dye tracing. *J Comp Neurol* 479:309-327. **NeuroVue Maroon** (previously PTIR271), **NeuroVue Red** (previously PTIR278)
- Hellard D, Brosenitsch T, Fritsch B, Katz DM. 2004. Cranial sensory neuron development in the absence of brain-derived neurotrophic factor in BDNF/Bax double null mice. *Dev Biol* 275:34-43. **NeuroVue Maroon** (previously PTIR271))
- Hsieh CY, Cramer KS. 2006. Deafferentation Induces Novel Axonal Projections in the Auditory Brainstem After Hearing Onset. *J Comp Neurol* 497: 589-599 **NeuroVue Red** was used for all figures except Figure 2D, for which both NeuroVue Red and Dil were used, and Figure 5A, for which Dil was used (K. Cramer, personal communication).
- Hsieh CY, Hong CT, Cramer KS. 2007. Deletion of EphA4 Enhances Deafferentation-Induced Ipsilateral Sprouting in Auditory Brainstem Projections. *J Comp Neurol* 504: 508-518. **NeuroVue Red**
- Jensen-Smith H, Gray B, Muirhead K, Ohlsson-Wilhelm B, Fritsch B. 2007. Long distance three-color neuronal tracing in fixed tissue using NeuroVue® dyes. **Immunol. Invest.**, 36, No 5-6, 763. **NeuroVues Jade, Red and Maroon**
- Maklad A, Fritsch B, Hansen LA. 2004. Innervation of the maxillary vibrissae in mice as revealed by anterograde and retrograde tract tracing. *Cell Tissue Res* 315:167-180. **NeuroVue Maroon** (previously PTIR271) was used for Figure 6 (B. Fritsch, personal communication)
- Matei V, Pauley S, Kaing S, Rowitch D, Beisel KW, Morris K, Jones K, Fritsch B. 2005. Smaller epithelia in Neurog1 null mice are related to earlier hair cell terminal mitosis. *Dev Dyn*, 234:633-650. **NeuroVue Maroon** (previously PTIR271) was used for Figure 8 (B. Fritsch, personal communication)
- Matei VA, Feng F, Pauley S, Beisel K, Nichols MG, Fritsch B. 2006 Near infrared laser illumination transforms the fluorescence absorbing X-Gal reaction product BCI into a transparent, yet brightly fluorescent substance. *Brain Res Bull.* 70:33-43. **NeuroVue Maroon** (previously PTIR271; see Fritsch, *Dev. Dyn.* 2005, above, for methods)
- Morris JK, Maklad A, Hansen LA, Feng F, Sorensen C, Lee KF, Macklin WB, Fritsch B. 2006. A disorganized innervation of the inner ear persists in the absence of ErbB2. *Brain Res.* 1091: 186-199 **NeuroVue Maroon, NeuroVue Red**
- Pauley S, Wright TJ, Pirvola U, Ornitz D, Beisel K, Fritsch B. 2003. Expression and function of FGF10 in mammalian inner ear development. *Dev Dyn* 227:203-215. **NeuroVue Maroon** (formerly PTIR271) rather than Dil was used in Figure 7 (B. Fritsch, personal communication)
- Pauley, S., Matei, V., Beisel, K.W., Fritsch, B. 2005. Wiring the Ear to the Brain: the Molecular Basis of Neurosensory Development, Differentiation, and Survival. In: *Development of the ear.* Eds. M. Kelley and D. Wu, SHAR, Springer Verlag, 85-120. **NeuroVue Maroon** (previously PTIR271) was used in Figure 4; see Tessarollo, *J Neurosci* 2004, below, for methods (B. Fritsch, personal communication)
- Tessarollo L, Coppola V, Fritsch B. 2004. NTF3 replacement with brain-derived neurotrophic factor redirects vestibular nerve fibers to the cochlea. *J Neurosci* 24:2575-2584. **NeuroVue Maroon** (previously PTIR271), **NeuroVue Red** (previously PTIR278)
- Zou D, Silvius D, Fritsch B, Xu PX. 2004. Eya1 and Six1 are essential for early steps of sensory neurogenesis in mammalian cranial placodes. *Development* 131:5561-5572. **NeuroVue Maroon** (previously PTIR271) and **NeuroVue Red** (previously PTIR278) were used for Figure 6, panels G-R (B. Fritsch, personal communication).

NeuroVue® is a trademark of PTI Research, Inc. used under license. NeuroVue® products are sold under sublicense from PTI Research, Inc.